

ADVAPIX TPX3 Version 1.0 - Datasheet

Model No.: APXMD3-Xxx170704












Illustration of single particle sensitivity of Timepix3 device. The tracks of different particles of radiation background (mostly muons and few protons) were recorded in 5 minutes on board of airplane. No noise (clean zero) is seen in dark regions. Inset shows the time profile along one muon track. The **ADVAPIX TPX3** modules were designed with special emphasis to performance and versatility which is often required in a scientific experimental work. They contain latest CERN detector Timepix3 for particle tracking and imaging. The fast modules with Si or CdTe pixel detectors Timepix3 can be used in different configurations such as stack of several layers or tiling to cover larger area. Each module contains single Timepix3 device with ultra-fast sparse data readout to acquire up to 40 Mhits per second. A separate USB 3.0 channel for each device assures fast read-out of the whole modular system. The sensor type and thickness is of customer's choice.

ΛΟΥΑΓΑΜ

Imaging the Unseen

The **ADVAPIX TPX3** can be used in a variety of applications such as energy resolved radiography (X-rays, neutrons, ions), particle tracking, time-of-flight imaging, multilayer Compton camera and many other. The sensors can be adapted for neutron imaging by deposition of converter layers¹. Recording individual hits together with advanced data processing allows increasing the spatial resolution in some applications to units of microns or even sub-micrometric level (for ions).

Compared to the older Timepix chip the Timepix3 detector offers 6 times better time resolution, 2 times better energy resolution, 2 times lower minimum energy threshold, zero dead time and 10x faster data transfer.

Main Features:

- Readout chip type..... Timepix3
- Pixel size55 x 55 μm
- Sensor resolution......256 x 256 pixels
- Time resolution.....1,6 ns
- Power..... External or via second USB 3.0
- Interface.....USB 3.0(Super-Speed)
- Maximum readout speed......40 million pixels / s
- Dimensions125 x 79 x 25.5 mm
- Weight 503 g

¹ Convertors based on 6LiF or 10B4C for slow neutrons (efficiency up to 4%) or PE for fast neutrons.





Device Parameters

Operating Conditions

Symbol	Parameter	Min	Тур	Max	Units	Comment
ТА	Temperature Range	0	50	70	°C	
Φ	Humidity			60	%	Not condensing

Electrical Specification

 $T_A = 25^{\circ}C$, USB voltage $V_{CC} = 4.8V$

Symbol	Parameter	Min	Тур	Max	Units	Comment
Vcc	Supply Voltage	4.0	5.0	5.5	V	
Icc	Supply Current					
ICC1	Chip disabled		800	1500	mA	
P 1	Power Dissipation			2.5	W	
I/O Conn. Inpu	ut CMOS (pin 4,6,8,10)					
VINL	Voltage Low	0		1.15	V	
VINH	Voltage High	2.15		3.3	V	
I/O Conn. Input LVDS (pin 3,5,7,9)						
VIN	Voltage Range	0		2.5	V	
VINDIFF	Differential Voltage	250		600	mA	
Bias Voltage S	Source for Sensor Diode					
VBIAS	Bias Voltage	0		±450	V	Polarity is sensor dependent


Performance Characteristics for Timepix3

Symbol	Parameter	Min	Тур	Max	Units	Comment
f	Hit-rate			40	MPixels/s	with USB 3.0 cable
	Data rate			2.4	Gbit/s	with USB 3.0 cable
TREAD	Frame Readout Time		33		ms	with USB 3.0 cable
dT	Time resolution	1.56			ns	
FREAD	Read-out frequency		320		MHz	1⁄2 of maximum ROC freq




Sensor parameters $T_A = 25^{\circ}C$

Symbol	Parameter		Si			CdTe	Units	Comment
	Thickness	100	300	500	675	1000	μm	
σ	Energy resolution of energy discrimination threshold (σ @ 23 keV)	0.5		1.1	keV			
σ	Energy resolution of energy discrimination threshold (σ @ 60 keV)	0.6			1.5	keV		
σ	Energy resolution in full spectral mode (σ @ 23 keV)	0.7			3.0	keV		
σ	Energy resolution in full spectral mode (σ @ 60 keV)		1	.0		3.6	keV	
	Typical detectable energy range for X-rays ²		1	.0		3 to 600	keV	See chart below
	Pixel size		5	5		55	μm	



² The maximum detectable energy is limited only by sensitivity of the selected sensor for the given radiation type. The maximum measured energy can reach several MeV in case of heavy charged particles



+5VDC connector

Main power supply (via standard 2.1 mm connector). Connect after plugging USB connector.

USB connector

USB type micro B, Standard USB 3.0 Super-Speed.

I/O Digital connector

Connector is not used in this version of the product.

1	GND		2	Max	
3	Reserved	LVDSoP (2.5V)	4	Reserved	CMOS o-5V
5	Reserved	LVDSoN (2.5V)	6	Reserved	CMOS 0-5V
7	Reserved	LVDS1P (2.5V)	8	Reserved	CMOS o-5V
9	Reserved	LVDS1N (2.5V)	10	Reserved	CMOS o-5V




Datasheet | Mechanical Dimensions



Mechanical Dimensions

Without protecting cover



All dimensions are in mm.

Extreme care must be taken when removing protecting cover and handling the **ADVAPIXTPX3** without the protecting cover. Warranty does not apply to mechanical damage of the sensor and wirebonds.








Datasheet | Configuration Examples



With protecting cover







All dimensions are in mm.



8

Ó

6

70





Model number codes

Device name: APX – AdvaPIX Device modification:		
	-	
MD3 – Timepix3 module		
Sensor type: P – Planar silicon E – Edgless silicon C – CdTe		
Sensor thickness: 1 – 100 μm 3 – 300 μm 5 – 500 μm		

Release history

Date	Changes
17/11/02	Model number codes added, datasheet version

Warning

Instructions for safe use. To avoid malfunction or damage to your **ADVAPIX** TPX₃ please observe the following:

• Do not expose to water, moisture.

• Do not disassemble. Wire-bonding connection may be irreversibly damaged.

• Do not insert any object into the sensor window.

Copyright

ADVACAM s.r.o.

U Pergamenky

1145/12,170 00

Czech Republic

Tel: +420 603 444 112 Email: info@advacam.cz

www.advacam.com

