

# WidePIX

CHROMATIC 5  
CHROMATIC 10

former name: WIDEPIX L<sub>2(1)x5</sub> - MPX3  
Model No.: WxAM3x-Xxx210721



Datasheet

## General features




Illustration of multichannel “color” radiographs where different materials are identified and imaged in different colors

The large area imaging detectors **WidePIX CHROMATIC 10** (Double Row) with resolution of 512 x 1280 pixels and **WidePIX CHROMATIC 5** (Single Row) with resolution of 256 x 1280 pixels are composed of Medipix3 hybrid detector electronics tiles. Each tile (256 x 256 pixels) is attached to a silicon or CdTe sensor. Therefore, the whole area of the **WidePIX CHROMATIC** device is fully sensitive and there are no gaps between sensor tiles. Each pixel has two integrated 12-bit digital counters and two energy discrimination thresholds. The counters store number of registered particles, e.g. X-ray photons, with energy above the appropriate threshold. Both counters can be joined to a single 24-bit counter providing enhanced dynamic range. The particle counting principle eliminates any additional noise generated by the sensor or electronic readout. It allows acquiring X-ray images with very high contrast and wide dynamic range. Therefore, even low contrast structures such as plastic or soft tissue are easily detectable in X-ray images.

Both devices are suitable for CT scanners, which can take advantage of large sensitive area without any gaps. The **WidePIX CHROMATIC 5** (Single Row) variant, moreover, supports a hardware-based Time-Delayed-Integration mode for online (continuous) scanning applications.

The energy discrimination thresholds of Medipix3 technology allow spectral X-ray imaging. Different materials in an inspected sample could be then identified based on their spectral X-ray attenuation properties. Energy spectra could be measured typically from 5 keV upwards.

The Charge Summing Mode implemented in the pixel electronics provides hardware-based correction of signal cross talk between pixels. This further considerably improves the detector’s spectral response and therefore also quality of spectra measured in individual pixels.

The camera is connected to a computer via an ethernet cable.

## Main Features

- Readout chip type..... Medipix3
- Pixel size<sup>1</sup>..... 55 x 55  $\mu\text{m}^2$
- Sensor resolution..... 512 (256) x 1280 pixels
- Dynamic range in one frame<sup>2</sup>..... 12-bit / 24-bit
- Dark current..... none
- Interface ..... 1x RJ45 1Gb/s ethernet
- Maximum frame rate<sup>2</sup>..... up to 80 (170) fps
- Dimensions ..... 170 x 140 x 42 mm
- Weight..... 2000 g

<sup>1</sup> 55 x 110  $\mu\text{m}^2$  at the edges and 110 x 110  $\mu\text{m}^2$  at the corners

<sup>2</sup> Depends on operation mode.

## Device parameters

### Operation conditions

Symbol	Parameter	Value	Units	Comment
T <sub>a</sub>	Operating ambient temperature range <sup>1</sup>	0-40	°C	
Φ	Humidity	< 60	%	Not condensing
IP	IP rating	IP50		

<sup>1</sup> With temperature stabilization – see the paragraph below.

### Water cooling interface

Temperature stabilization of the device required when in operation. **WidePIX CHROMATIC** uses water connectors that allow for quick disconnection/reconnection. Mating connector is included as standard accessories and must be attached to 4x6 mm plastic hose.



Temperature of the cooling water must be within range  $21 \pm 4$  °C.  
Max. pressure in the water-cooling system: 1,2 bar.  
The device will automatically shut down after chip or CPU temperature exceeds 55 °C.  
Intended for dust free indoor use.

### Electrical specification

T<sub>a</sub> = 25 °C, V<sub>CC</sub> = 12 V

Symbol	Parameter	WidePIX CHROMATIC 5 (Single Row)	WidePIX CHROMATIC 10 (Double Row)	Units	Comment
V <sub>CC</sub>	Supply Voltage	12 ±10 %		V	
I <sub>CC</sub>	Supply Current (V <sub>CC</sub> = 12V)	0,7/1,4	1,6/3,2	A	Typ/Max
P	Power dissipation	9/18	19/38	W	Typ/Max
A	Sensor Area	70,5 x 14,1	70,5 x 28,2	mm	
	Detector Resolution	1280 x 256	1280 x 512	Pixels	
f	Frame Rate <sup>1</sup>	170	90	fps	
T <sub>READ</sub>	Readout Time <sup>2</sup>	6	13	ms	
m	Weight	1900	2000	g	

<sup>1</sup> Operating parameters: Shutter time = 1 ms, Mode = CSM or SPM-1Ch 12bit resolution.

<sup>2</sup> During Readout time (or Dead time), no charge is collected from the sensor.

## Threshold setting

Typical values for 300  $\mu\text{m}$  Silicon sensor,  $T_a = 20^\circ\text{C}$ .

Range	Mode	Min. energy threshold [keV]
Super Narrow	SPM	-
	CSM	Not supported
Narrow <sup>1</sup>	SPM	5,0
	CSM	5,0
Broad	SPM	Not supported
	CSM	Not supported

Typical values for 1000  $\mu\text{m}$  CdTe sensor,  $T_a = 20^\circ\text{C}$

Range	Mode	Min. energy threshold [keV]
Super Narrow	SPM	-
	CSM	-
Narrow <sup>1</sup>	SPM	8,0
	CSM	8,0
Broad	SPM	-
	CSM	-

## Electrical specification


$T_{\text{dev}} = 22^\circ\text{C}$

Parameter	Si		CdTe	Units	Comment
Thickness	300	500	1000	$\mu\text{m}$	
Bias Voltage	150	150	- 450	V	Max
Minimum energy threshold	5,0		8,0	keV	Typical, $T_A=22^\circ\text{C}$
Typical detectable energy range for X-rays <sup>2</sup>	up to 60		up to 600	keV	See chart below
Pixel size <sup>3</sup>	55 x 55		55 x 55	$\mu\text{m}^2$	

<sup>1</sup> By default the detector will be calibrated for the Narrow Gain Mode. Additional gain modes can be added upon request. Broad range and combination of Super Narrow mode with CSM is not available with silicon detectors.

<sup>2</sup> To get true detector response, detectable energy and quantum efficiency of sensor chip must be combined with energy range of readout chip.

<sup>3</sup> 55 x 110  $\mu\text{m}^2$  at the edges and 110 x 110  $\mu\text{m}^2$  at the corners.




## Basic principles, measurement types and operational modes

The ionizing radiation particle interacts with the sensor material creating an electric charge. This charge is collected by electric field and brought to pixel preamplifier where it is amplified and shaped forming triangular voltage pulse. The amplitude and duration of this pulse is proportional to energy deposited by particle within the pixel. The situation when the voltage pulse amplitude in particular pixel exceeds preselected threshold value is called “event” or “hit”.

Each pixel contains two digital counters (12 and 12 bits). These counters are used differently according to measurement type and mode. List of operational modes and their description is provided in the table below.

Type	Mode	Bit depth	Description
Frame (reading all pixels)	SPM-1CH	12/24 bit/frame	Single Pixel Mode using one counter: Every pixel works independently of its neighbors. One energy threshold (energy channel) is available. 1 output image: Number of events per pixel
	SPM-2CH	12 bit/frame	Single Pixel Mode using both counters: Every pixel works independently of its neighbors. Two energy thresholds (energy channels) are available. 2 output images: Number of events per pixel
	CSM	12/24 bit/frame	Charge Summing Mode: The charge from 4 adjacent pixels is summed and is assigned to the pixel with the largest charge deposition. The event is counted only if the sum of signals exceeds the second energy threshold. 1 output image: Number of events per pixel

Device description



Ethernet connector

1 x RJ45 1 Gbit/s ethernet connector

+12 V DC connector

Main power supply (via standard M8 connector with 3 female contacts). Connect after plugging ethernet cable. Connect the power adapter to the device first, then to the wall socket.

Synchronization interface (optional)

Two 4-pin M8 connectors (female for outputs and male for input) serve as synchronization interface, allowing to synchronize **WidePIX CHROMATIC** detector with external processes. Four signals are available:


- **Ready in** – measurement is not possible, when signal at logical zero
- **Trigger in** – logical zero starts shutter (measurement)
- **Ready out** – logical one if device is ready to for new shutter
- **Trigger out** – mirrors shutter (logical zero when shutter is active)

All signals are TTL compatible and 5 V tolerant. For a detailed description see **Synchronization Guide for WidePIX Chromatic** in the WidePIX Chromatic Manual.

Sync. Outputs (M8-4Female)		Sync. Inputs (M8-4Male)	
Pin	Signal	Pin	Signal
1	Gnd	1	Gnd
2	Trigger Out	2	Trigger In
3	Ready Out	3	Ready In
4	Reserved	4	Reserved

## Mechanical dimensions


The following drawings illustrate the standard connection setup of **WidePIX CHROMATIC 10** (Double Row) and **WidePIX CHROMATIC 5** (Single Row).




All dimensions are in mm.

## Sensitive area


Sensitive area for models **WUAM3B**-XxxYYMMDD (single upper row)



Sensitive area for models **WLAM3B**-XxxYYMMDD (single lower row)



Sensitive area for models **W2AM3B**-XxxYYMMDD (two rows)





## Instructions for safe use



**Do not touch sensor surface!**

To avoid malfunction or damage to your **WidePIX CHROMATIC** please obey the following:

- Do not expose to water or moisture **WidePIX CHROMATIC** is dust protected only.
- Do not open **WidePIX CHROMATIC** case. Detector wire-bonding connections may be irreversibly damaged.
- Do not operate detector when not properly water cooled. Otherwise, detector temperature may rise above the specified range. Recommended temperature is 22 °C.
- The protection provided by this product may be impaired if it is used in a manner not described in this document.

## Disposal



Do not dispose of these instruments as unsorted municipal waste. Please use separate collection facilities or contact the supplier from whom the instrument was purchased. Ensure discarded electrical waste is properly recycled to reduce environmental impact.

## Release history

Date (YY/MM/DD)	Changes	Changed by
19/07/28	Preliminary version	
21/06/26	ETH version	
23/03/14	New drawings and corrected versions	
23/08/08	Supply Voltage changed from 24 V to 12 V	
23/09/05	Default gain mode added	
24/02/15	Datasheet revision	J. Baborák
24/06/28	Water cooling details added	J. Baborák
24/07/02	New graphic style of the document	P. Bloudek
24/07/23	Minor format changes	J. Baborák
24/10/23	IP reset pointer added	J. Baborák
26/01/08	New graphic style of the document, rebranding and revision	P. Bloudek, S. Valtera

### **ADVACAM s.r.o.**

U Pergamenky 12  
Prague 170 00  
Czech Republic

Tel.: +420 608 605 533  
Email: [sales@advacam.com](mailto:sales@advacam.com)  
[www.advacam.com](http://www.advacam.com)

